



## Aeon Medical UK

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### Abdominal Drainage Kit

Specially designed, competitively priced abdominal kit offered by us have the following characteristics:

- Soft abdominal drainage catheter with atraumatic rounded open distal end.
- Radio opaque line on catheter
- Collection bag capacity 2000 ml
- E.T.O. Sterile & Pyrogen Free.

**Packing** : 10 pcs in each box.

**Size** : 16,20,24,28,32,36 FG



## Aeon Medical UK

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*