



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Anaesthesia Face Mask (Black Rubber)

Flexible 22 F ISO connection, transparent cushioned and crystal clear Anaesthesia Face Mask (Black Rubber)

- Air filled cushion to provide safety and comfortable use.
- Suitable use in respiratory care of anaesthesia and resuscitator.
- Body with antistatic property.
- Excellent grip.

Material used : Black Rubber

Packing : 10 pcs in each box

Size : All pediatric and adult 0, 1, 2, 3, 4, 5



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com