



## **Aeon Medical UK**

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### **Anaesthesia Face Mask (Silicon)**

Specially designed anaesthesia face mask (silicon) have a transparent cushioned, crystal clear, flexible 22 F ISO connection.

- Air filled cushion along with silicon is used to provide safety and comfortable use
- Suitable for use in respiratory care of resuscitator and anaesthesia.
- Body with anti static property.
- Good grip.

**Material used** : Silicon

**Packing** : 10 pcs each box.

**Size** : All pediatric and adult 0, 1, 2 , 3, 4, 5



## **Aeon Medical UK**

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*