



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Anaesthesia Face Mask

We offer crystal clear, transparent cushioned, flexible 22 F ISO connection anaesthesia face mask. These masks have the following features:

- Air filled cushion used to offer more comfort and safety.
- Suitable for use in respiratory care of anaesthesia and resuscitator.
- Body with anti static property.
- Comes with effective grip.

Packing : 10 pcs per box.

Size : All pediatric and adult 0, 1, 2 , 3, 4, 5



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com