



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Arterial Line Filter

The arterial line filter manufactured by us is completely based on new construction principles that help in increasing the efficiency of the filter and maximizing safety to the patient. The Arterial filter has the various features that gives it an edge over the products available in market are:

- Easy handling,
- Effective elimination of gas emboli, fat emboli
- Aggregates composed of Platelets, Red Blood Cells and other debris.
- Filtration membrane 40 micron.
- Priming Volume: 165 ml.
- Maximum Flow 7 L/min.
- E.T.O. Sterile & Pyrogen Free.

Packing : Box of 5 pcs.



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com