



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



AV Tubing

- Raw material: Non-toxic medical grade P.V.C. tubing.
- Vented protective caps are provided at tubing ends.
- Available in different lengths as per customer's requirement.
- E.T.O. Sterile & Pyrogen Free.

Packing : 10 pcs per box

Size : 1/2" X 3/32", 3/8"X 3/32", 1/4" X 1/16" Other size



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com