



## Aeon Medical UK

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### Belflon (I.V. Cannula)

The Belflon I. V. Cannula has been exclusively designed keeping in mind the needs of the clients. They are available with / without wings and an injection part to assist infusion therapy amongst neonates. 6% luer taper which are at par with any other international product available in market. It ensures pain free acceptable teflon made radiopaque shaft with metallic stillete.

- Non-toxic
- Injection part with one-way valve has a flange attached to the catheter / shaft
- Standard size hub attached to the distal end for intravenous line attachment.
- Injection part with color-coded cap.
- Color coding for different sizes.
- Transparent blister packing with tyvec paper.
- E.T.O. Sterile
- Pyrogen free.

**Packing** : 100 pcs per box.

**Sizes** : 14, 16, 18, 20, 22, 24g



## Aeon Medical UK

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : info@aeonmedicaluk.com*