



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



BIPAP & CPAP Circuit

Bipap circuits are available with 1.6/1.8/2.0 meter long corrugated flexible tubing with attached 22m/22f & 22m/22m adaptor both ends.

CPAP circuit are available with 2.5 meter long corrugated flexible tubing with attached 22m/22f & 22m/22m adaptor both ends & 22m/15f patient t-piece adaptor attached at 0.5m from one side.

They are individually packed in peelable pouch and have E.T.O. Sterile and Pyrogen Free.

Packing : 10 pcs. Per box



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com