




Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Cardioplegia Adapter - 4 Way

The Multiple Perfusen set supplied by us helps in simultaneous perfusion of the aortic root and up to three or more vein grafts. The set features an inlet port with female luer and clamp attached to an adapter which splits into four legs. The legs have locking male luers locks and color coded clamps for easy identification.

Packing : 10 pcs. per box



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com