



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Chest Drainage Catheter

- Manufactured from non-toxic medical grade P.V.C. tubing.
- To maximize drainage from the chest, distal end is open with super smooth large eyes.
- To provide excellent forceps grip for traumatic passage through tissue, specially designed tapered tongue fixed at proximal end.
- Catheters are marked at 2 cm from the last eye to ascertain the depth of placement.
- Radio - opaque individually double packed in peelable pouch.
- Angled catheters are also available. E.T.O. Sterile and Pyrogen Free.

Packing : Box of 25 Pcs.

Size : 16,20,24,28,32,36 FG



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com