



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Control Syringes

Single use control syringes offered by us are individually packed. They have a control ring for permanent binding to the syringe. The syringes offered by us have a polypropylene barrel, plunger rod and latex-free plunger tip.

Size : 10, 12, 20ml



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com