



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



CVP Manometer

- Suitable for continuous or intermittent monitoring of central venous pressure during infusion.
- Manometer tube with graduation from -4cm to +34cm and attached to three way connector.
- 125 cm long extension tube is fitted with female luer lock at one end and male luer lock on the other end for safe connection to C.V.P. Catheter.
- Latex rubber flash ball is provided for extra medication.
- Sliding indicator provided to record previous reading.
- Moulded clamps with latex rubber strings are provided for easy fixing on I.V. Stand.
- Individually packed in peelable pouch.
- E . T . O . Sterile and pyrogen free

Packing : 40 pieces of boxes



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com