



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Dialysis Catheter Kit

We offer three different types of Hemodialysis Catheters, namely, Straight Catheter, Bent Catheter and Curved Extension Catheter. These hemodialysis catheters offered by us are available in different packings offered by the clients. We offer these kits in special individual packing or in combinations in trays.

Hemodialysis Catheter Kit comprises of :

- Hemodialysis Catheter (this is Polyurethane Catheter, Radiopaque, with length makers, color coded safety clamps, luer connection hubs and protective caps)
- Scalpel No. 11
- Introducer Needle.
- Luer Lock syringe , 5MI
- Tissue dialator.
- Luer Lock Injection Caps.
- Guide Wire (J Tip) in dispenser, with protection cap for the guide wire tip.

Combinations for hemodialysis catheters



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com