



## Aeon Medical UK

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### Disposable Chest Drainage System (Adult)

- The surgical chest drainage bottle has 2.5 liter capacity making it suitable for cardiac and thoracic procedure.
- Easy to read graduation helps in determining the drain volume precisely.
- Clearly marked initial level ensure the under water seal.
- Write on facility on the side of scale helps monitoring of drain volume.
- Separate suction port is provided to connect with suction unit.
- Kink resistant large bore tubing facilitates unrestricted flow, unit is provided with hanger and floor stand.
- Complete unit is sterile and packed in a box. E.T.O. Sterile and Pyrogen Free.

**Packing** : 5 pcs per box.



## Aeon Medical UK

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*