



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Disposable Suction Bag

Disposable Suction Bag offered by us is first of its kind, integrated overflow protection / bacterial filter, odour protection due to laminated film and adaptable to suction units or any other pipeline wall suction system.

Capacity : 1 ltr and 2 ltr



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com