



## **Aeon Medical UK**

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### **Endotracheal Tube (Cuffed)**

Our Endotracheal Tube (Cuffed) is manufactured from siliconised non-toxic and non-irritant P.V.C. kink resistant thermo sensitive tube that ensures the tube patency for patient safety and abreast of it softens at body temperature to conform to the anatomy of the respiratory tract. The features are as follows:

- The murphy eye has a polished smoothness that ensures atraumatic intubation and extubation.
- During cuff inflation and deflation, the kink resistant inflation tube ensures patient safety.
- Between the tube and the trachea, the soft low pressure cuff provides an air tight fit.
- The X-ray opaque line and markings that are provided on the tube helps to know the exact location of the tube position.
- A standard 15mm connector is provided at the proximal end. E.T.O. Sterile & Pyrogen Free.

**Packing** : Available in box of 10pcs.

**Size** : 6.0/6.5/7.0/7.5/8.0/8.5/9.0/9.5/10 MM



## **Aeon Medical UK**

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*