



## **Aeon Medical UK**

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### **Endotracheal Tube (Plain)**

Endotracheal Tube made from special non-toxic medical grade material

- It is clear and made of thermo sensitive PVC material.
- It has radio opaque line to guide for the exact location of tube.
- All the Endotracheal Tube is fitted into it with 15mm standard connector..
- The smooth tip helps atraumatic intubation.
- It is E.T.O Sterile and Pyrogen Free.

**Packing** : Available in box of 25pcs.

**Size** : 2.5/3.0/3.5/4.0/4.5/5.0/5.5/6.0/6.5/7.0 MM



## **Aeon Medical UK**

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*