



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Heat & Moisture Exchange Filter

- Tidal volume range (ml): 70-600
- Humidification efficiency (mg/l)
 - at Tidal Volume = 100 ml:33.0
 - at Tidal Volume = 250 ml:31.0
 - at Tidal Volume = 500 ml:27.5
 - at Tidal Volume = 1000 ml
- Bacterial Filtration (%):99.98
- Weight (g):16.5
- Viral Filtration (%) :99.80
- Compressible volume (Dead space) (ml):22
- Resistance (cm H₂O) at flow:
 - 0.15 l/s: 0.6
 - 0.25 l/s: 1.0
 - 0.50 l/s: 2.1
 - 1.0 l/s
 - 1.5 l/s

Packing : Box of 100 pcs.

Size : Adult , Paediatric and Infant



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com