



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Hemodialysis Blood Tubing Set

The specially designed blood tubing set are used for connecting external system which extracts patient's blood to the dialyzer and then returns patient's blood from the dialyzer. These can be customized as per the requirements of the client's.

- It consists of two parts - arterial and venous, which are used during dialysis and attached with fistula and dialyzer.
- The unique chambers reduce foaming, increase air removal and do not trap epo compliant material of transducer protectors. This helps in ensuring secure machine fit and less incidence of wet out.
- Color coded ergonomic clamps and dialyzer connectors with tethered caps that seal tightly for disposal.
- Kink resistant blood tubing is easy to handle and ensures fast priming.
- E.T.O. Sterile & Pyrogen Free.

Packing : 10 pcs in each box.



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com