



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



I.V. Line Filter

- Self & fast priming of adult I.V. Set to minimize nursing time.
- Able to withstand a minimum of 45 psi (3.1 bar) pressure.
- Automatic air venting in any position.
- Loss of flow due to air locks.
- Bacterial retention throughout 96 hour flow.
- Filter media is 0.2 micron supor.
- E.T.O. Sterile and Pyrogen Free.

Packing : 25 pcs per box



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com