



## **Aeon Medical UK**

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### **Inter Coastal Drainage System**

---

- Helps in collection of drainage fluid from thoracic and abdominal cavity
- Exclusively designed moulded handle for easy carrying and hanging for bag
- Volume of PVC drainage bag: 1000 ML capacity
- Soft and kink resistant non-toxic medical grade PVC tubing
- E.T.O. Sterile and Pyrogen Free.

**Packing** : 10 pcs. per box



## **Aeon Medical UK**

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*