



## **Aeon Medical UK**

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### **Levins Tube**

- Suitable for gastro intestinal feeding and aspiration.
- Distal end is coned with four lateral eyes.
- Radio opaque line throughout the length.
- Low friction medical grade P.V.C. reduces the incidence of nasal discharge.
- E.T.O. Sterile and Pyrogen Free.
- Proximal end is provided with universal funnel.
- Marked at 45, 55, 65 & 75 cm from the distal end.

**Packing** : Box of 25 pcs.

**Size** : FG 8/10/12/14/16/18/20



## **Aeon Medical UK**

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*