



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Measured Volume Burette Set

- Large calibrate burette chamber.
- E.T.O. Sterile and Pyrogen Free.
- Injection port for supplementary medication.
- Automatic shut off valve, Micro Dropper, Flexible drip chamber
- Luer lock connector also available.
- Bacteria proof venting.
- Sharp piercing spike provided with build-in air inlet containing bacteria barrier filter.
- Roller type flow controller.
- Y injection sites for extra medication.

Packing : 25 pcs per box

Size : Chamber Capacity: 110ml / 150ml



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com