



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Mini Drain - O - Pack

Mini Drain - O – Pack are very suitable for laproscopic surgery. Jackson plate type 7mm Flat drain is manufactured from extra soft PVC, X-ray opaque line provided throughout the length of catheter. E.T.O. Sterile and Pyrogen Free.

Packing : Box of 20 pcs.

Size : 6 FG/8 FG



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com