



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Scalp Vein Set

- Butterfly set for long term infusion
- Short bevelled siliconised needle facilitates atraumatic cannulation.
- Rigid luer lock fitting also available on specific demand.
- Thin wall needle provides better flow rate per gauge, as in the same cannula more circumference is available for better flow of fluid.
- Butterfly shaped wings facilitate easy handling and attachment with the skin.
- The proximal end of the set is provided with flexible female luer fitting.
- Renders rapid venous access
- Greater patient comfort during infusion.
- Butterflies are color coded for instant identification of needle size.
- Butterfly is connected to soft non toxic, non-irritant, medical grade tube, which does not kink or coil set.
- E.T.O. Sterile & Pyrogen Free.

Packing : Box of 100 pcs.

Size : 18,19,20,21,22,23,24,25



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com