




Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Silicon Venteuse Cup

Silicon venteuse cup helps in safe birth assistant which minimizes risk of trauma and cephalohematome of fetal cranium. It is available with sterilizable cup of 60 mm & 70 mm each. It is connectible to any suction unit with gauge & release valve.



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com