



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Spinal Needle

Our superior quality spindle needle is offered with plastic hub and quincke point. The spindle needles that we provide are available with unique translucent window hub, which allows for easy visualization. The important features of our spindle needle are as cited below

- It has specialized fitted stylet to reduce tissues coring
- It is E. T. O. sterile and pyrogen free

Length : 1.5 & 3.5

Packing : Available in box of 25 pieces



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com