



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Suprubic Balloon Catheter with Trocar Cannula

Specially designed for cystostomy and supra – pubic drainage and made of non-toxic, non irritant medical grade P.V.C.

- For maximum drainage: one lateral eye at distal end
- Symmetrical inflated balloon provided for self-retention
- Balloon capacity 10 ml
- Radio opaque line.
- E.T.O. Sterile & Pyrogen Free.

Packing : 10 pcs per box.

Size : 16/18/20/22 FG



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com