



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Syringes with Needle

- Prominent graduation ensures accuracy of dosage and avoids wastage of the medicine.
- Clear transparent chamber comparable to glass syringe facilitates visual inspection.
- Thermoplastic elastomer gasket is inert for minimum friction during movement and prevents leakage and back flow.
- Sharp triple faceted bevelled tip needle for minimum trauma during penetration.
- Ribbon packed in tough laminated film
- Highly resistant to puncture and tearing
- E.T.O. Sterile & Pyrogen Free. Ready for use.

Size : 1 ml, 2 ml, 3 ml, 4 ml, 5 ml, 10 ml, 20 ml, 50 ml



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com