



## **Aeon Medical UK**

*suppliers of surgical disposables*

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*



### **Tracheostomy Tube (Cuffed)**

- The non-return valve along with indicator facilitates the cuff inflation and deflation.
- The terminal end of the tube is suitable for 15mm connector.
- In peelable pack sterile is supplied with a lid for convenience.
- The thin smooth cuff of the tube inflates symmetrically and eliminates any risk of herniation.
- It also ensures the accurate positioning of the tube in the trachea.
- It is E.T.O. Sterile & Pyrogen Free.

**Packing :** Available in box of 25 pcs  
**Size :** 5.0/6.0/7.0/7.5/8.0/9.0/10.0 MM



## **Aeon Medical UK**

*ISO 9001:2000, CE Certified, WHO GMP Compliant Company*

*Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.*

*Tel. : +441277201166 Fax : +441277200347*

*Email : [info@aeonmedicaluk.com](mailto:info@aeonmedicaluk.com)*