



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Umbilical Catheter

- Umbilical Catheter are designed for intermittent or continuous access to the umbilical artery or vein of newly born or premature infant.
- Made from non-toxic soft medical grade P.V.C. that yields easily to tissue contours at body temperature.
- Individually straight packed in peel-off pouch packing to ensure aseptic handling.
- E.T.O. Sterile & Pyrogen Free.

Packing : 25 pcs per .

Size : 3.5 / 5 / 6 / 8 FG



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com