



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Umbilical Cord Clamp

Fabricated for intermittent or continuous access to the vein or umbilical artery of newly born or premature infant. Manufactured from non-toxic soft medical grade P.V.C. that yields easily to tissue contours at body temperature. Each and every clamp is straight packed in peel-off pouch packing to ensure aseptic handling. E.T.O. Sterile & Pyrogen Free.

Packing : 25 pcs in each box

Size Available : 3.5 / 5 / 6 / 8 FG



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com