



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Urine Bag With Measuring Cup

- Specially designed for accurate measurement of hourly urine output.
- Accurate measurement.
- Complete close circuit eliminates the risk of contamination.
- 2000 ml urine bag attached directly with measured volume meter ensures auto over flow
- Convenience to empty the measured volume meter in bag.
- E.T.O. Sterile and Pyrogen Free.

Packing : 10 pcs. Per box



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com