



Aeon Medical UK

suppliers of surgical disposables

ISO 9001:2000, CE Certified, WHO GMP Compliant Company



Ventilator Circuit

- Medical Grade Corrugated PVC ventilator tube circuit can be ported for temperature and capnography, circuit is provided with a recessed standard 22mm Male / Male & 22 mm Male / Female adapter for flexibility throughout the critical care management.
- The unique water trapper offers excellent drainage and ideal positioning to maximize moisture collection.
- It can be assembled according to the customer requirement.
- E.T.O. Sterile and Pyrogen Free.
- Self expandable and non-expandable (without water trapper / with one water trapper / with two water trapper)

Packing : Box of 20 pcs.

Size : Adult & Paediatric

Color : Blue and Clear



Aeon Medical UK

ISO 9001:2000, CE Certified, WHO GMP Compliant Company

Aeon Medical UK Ltd., Linden Lea, Hillwood Grove, Hutton Mount, Essex CM 13 2 PD, UK.

Tel. : +441277201166 Fax : +441277200347

Email : info@aeonmedicaluk.com